



“The Against the Gravity”

# N-COG Lift introduction and USER GUIDE

[www.nfinders.com](http://www.nfinders.com)

**N-FINDERS**

**N-COG**  
엔-코그 [Newton Cog]

엔스캐폴드  
**N-scaffold**

N-COG  
**Royal**

N-COG  
**Y-KO**

N-COG  
**Spiral**

+ 더마스프링  
**Derma Spring**

+ 더마스프링 아이  
**Derma Spring eye**

Copyright © 2012 by N-Finders Co., Ltd.

This document is solely for the use of client personnel.

No part of it may be circulated, quoted, or reproduced for distribution outside the client organization without prior written approval from N-finders Co., Ltd.


This document provides an outline of a presentation and is incomplete without the accompanying oral commentary and discussion.

# Agenda

1. What is Polydioxanone? (PDO)
2. Minimally invasive lifting with threads
3. Effects and Clinical importance
4. Why N-FINDERS?
5. Type of threads
6. Procedures with threads
7. Avoiding complications
8. Closing



# What is PDO?



Polydioxanone chemical structure

- Mainly use for biomedical applications
- Colorless, crystalline, a biodegradable synthetic polymers
- 100% biodegradable material which can be completely absorbed by the body and excreted through urine, exhaled as CO<sub>2</sub> entirely after 6-8 months.
- Excellent flexibility, safe treatment without acute inflammation reaction.
- Used in surgical sutures and orthopedic

# Minimal Invasive Lifting with Threads


- PDO thread lifting treatment tightens sagging skin tissue instantly
- The threads are inserted by needles, anchored, and pulled to lift the skin upwards and tightening the skin of the face

**Instant lifting effect**

+

**Long term effect : Neocollagenesis , Neovascularization**

- Procedure takes only 15 - 30 minutes





# Main Effects

- Instant Skin Lifting
- Collagen stimulation -> Neocollagenesis
- Neovascularization
- Increase of growth factors synthesis
- Cellulite decrease



**Skin elasticity**



**Skin Tightening, by contracting fat tissue**

# Clinical Importance

- Just local anesthesia: you can return to normal life immediately after the procedure.
- Minimal downtime
- No incisions or stitches are required
- No scars
- No secondary effects
- Suitable for everyone and applicable to the whole body

Why

# N-FINDERS

- Scientific background
- Clinical Background
- Educational system
- Several types of threads (and more coming!)
- **Hand Made (molded thread x)**

ISO 9001, ISO 13485, KGMP, CE, FSC, MFDS Certificate

Certificate of trademark registration/ ISO 9001:2015 / ISO 13485:2016 / EC Design Examination certificate/

Certificate of service mark registration / Certificate of patent



Certificate of Manufacturer for  
Mono/GSI company Prefix number  
license/ Certificate of  
GMP/Certificate of manufacturer for  
Cogl

Certificate of trademark  
registration/ ISO 9001:2015 /  
ISO 13485:2016 / EC Design  
Examination certificate/  
Certificate of service mark  
registration / Certificate of patent  
/



# Why N-FINDERS

## EVIDENCE BASED MEDICINE

N-Cog lifting's effect has tested and proven by medical experiments, and clinical data and scientific papers (SCI reports) published on international journals.




Various SCI reports



### Derma-Spring Study highlight Bulge formation after inserting spring type PDO thread



(L) Insertion of 3-5 spring type PDO thread shows augmentation effect.



(R) Inside the tissue, the spring shape remains almost as same as it was before.



### N-Cog Study highlight Bi-directional, Multi-directional Cogged Thread

#### Bi-Directional PDO thread insertion



Skin is folding in accordance with inserted threads

#### Multi-directional PDO thread insertion



Folding is not observed





# Why N-FINDERS

N-Finders has published several **PDO threads training books**, written by professors of human anatomy, plastic surgeons and dermatologist


A one of a kind complete and total guide series for perfect PDO threads lifting



# N-Cog PDO Threads

## ○ MONO THREAD

- ☐ RE:MONO
- ☐ SCREW
- ☐ DERMA SPRING
- ☐ DERMA SPRING EYE
- ☐ SCAFFOLD



## COG THREAD

- ☐ N-COG
- ☐ SPIRAL
- ☐ ROYAL
- ☐ Y-KO
- ☐ N-FLIP



## Micro Punching Thread

- ☐ N-FIX



# RE:Mono

- Skin tightening and cellulitis improvement
- Prevention of sagging
- Volumizing

RE:




Screw

SCREW  
THREAD



Derma-Spring

+ 더마스프링  
**Derma  
Spring**



Derma-Spring EYE

+ 더마스프링 아이 eye  
**Derma  
Spring**



# Cog threads:N-Cog






# Cog Threads: N-Cog

- RE: N-COG MULTI is a multiple bidirectional barbed PDO, specially designed to prevent skin dimpling, protrusion and breakage of the thread once it is inserted inside the skin tissue.

## N-Cog Multi



Needle (mm)		Thread(mm)	
Gauge	Length	USP	Length
23G	60	3-0	100
23G	90	3-0	130
Half-Cut cannula			
21G	60	2-0	100
21G	100	2-0	150

30 threads per box. (6 pouches with 5 units in individual sachet)

## TREATMENT PROCESS



## INSTRUMENT READY



- \*An array of surgical instruments for thread lifting /1 patient
- (A) 10cc syringe 1ea
- 30G 38mm long needle 5ea
- (C) Alcohol swab-many
- (D) Sodium chloride injection 20 ml 2ea
- (E) Epinephrine 1ml - , Lidocaine HCL Hydrate 2%, 20ml / Normal saline 100ml , sodium bicarbonate 5cc
- OR >Xylocaine 2% with adreanaline 2ea
- (F) Operation drape (60cm\*60cm) 1ea
- (G) Tower clip
- (H) IRIS operating scissors
- (I) IRIS forceps
- (J) Dressing forceps
- (K) Sterilized latex surgical Glove
- (L) Dressing package Tray
- (M) Emesis Basin
- (N) Tray cover
- \*Puncture needle 18G, 25G or needed one
- \*Dental Lidocaine –in case need
- \*Dental Lidocaine injector – in case need

## TREATMENT PROCESS

### Design

Draw a design line on the face or treatment part  
Considering patient's condition



## TREATMENT PROCESS

### ANESTHESIA - Anesthetic cream

9~11% lidocaine cream for 30~40 min  
Wipe the cream after





## TREATMENT PROCESS

### ANESTHESIA - Local anesthesia

Inject the anesthesia solution evenly in the area to be treated using a long needle. This is the basic anesthesia solution and it will be made according to user's preference.

Normal saline 50ml  
+ Lidocaine HCL Hydrate 2%, 20ml (anesthetics)  
+ sodium bicarbonate,(Na -bicarbonate) 5cc  
+ Epinephrine 1ml




## TREATMENT PROCESS

### ANESTHESIA- Nerve block


\* 30G needle with 10cc syringe with anesthetic solution (ex: lidocaine + Epi ,1% lidocaine 1~2cc injection )

\* Entry point :: local anesthesia with dental lidocaine




1. Midface : Infraorbital nerve - the foramen lies along a line parallel to the midline that passes through the pupil and corner of the mouth.


- Area : lateral part of nose, lower eyelid, upper lip, cheek



2. Lower face : Mental nerve  
- Area : chin, lower lip



3. Upper face : supraorbital nerve  
- Area : supraorbital, supratrochlear, upper eyelid, forehead



## POSSIBLE PRESCRIPTION

Prescriptions must be made according to the patient's circumstances at the judgment of the attending physician, and the responsibility for all prescriptions rests with the attending physician.

1. Aspirin : should not be terminated. However, bruising is more likely.
2. Patients at risk of Herpes labialis : advised to take prophylactic acyclovir tablets and should apply topical cream, especially if threads are to be inserted around the lips.
3. Arnica tablets : taken from about 5 days pre-treatment may be beneficial in reducing swelling after the treatment.
4. all blood 'thinning" supplements ( Evening Primrose Oil, Vitamin C and Vitamin E and certain herbs) : Stop taking



## GENERAL AFTERCARE ADVICE

Some pain/discomfort after treatment can be expected, Ibuprofen can be taken if required

Do not have any other aesthetic procedures (such as toxins, fillers, peels etc.) for 7 days

Avoid extreme temperatures, alcohol and anticoagulants for 7 days. Those taking Prescribed anticoagulants such as aspirin/warfarin must continue to do so but must be aware that bruising is more likely

Gently wash the treated area as required, but do not rub or massage for 2 weeks

It is important not to overstretch the treated area (i.e. open the mouth wide, avoid dental treatment about 2weeks if possible)

Do not have Radio-frequency, IPL, Laser or other heat treatments near the treated area for 2 weeks



Do not go to the sauna or do vigorous exercise on the day of the procedure

## POSSIBLE SIDE EFFECTS

1. **Common and mild side effects are possible from PDO Lift** : Swelling, Soreness, Bruising, Temporary skin dimpling, mild side effects could last 2~3 days.
2. **Patient may continue normal activities immediately.**
3. The side effects are mostly manageable – removal of threads may be required in cases such as migration and extrusion. There was no scarring after removal, and the problems resolved completely.
4. Serious but very rare : Infection, Granuloma, Protrusion, Facial nerve trauma, Asymmetry



# Procedures : N-Cog



## MID FACE 07-2

**Jowl : N-Fix, N-Cog Multi / N-Cog Spiral**


**Entrance point :** ● Temporal hair line

**End point :** Extension line from mouth corner to ear lobe


**Layer :** SMAS layer in temple and cheek



# Procedures : N-Cog




# Procedures : N-Cog






# Procedures : N-Cog






# Procedures : N-Cog



# Procedures : N-Cog



## LOWER FACE 13-1

### Marionette line : N-Cog Multi / N-Cog Spiral

**Entrance point :** ● 1 cm anterior to the tragus

**End point :** 1 cm medial to the marionette line

**Layer :** SMAS layer



N-Cog Multi 21GX100 mm



or




N-Cog Spiral 21GX100 mm

To improve marionette lines, mark the marionette line with a marker and insert 3-5 strands of N-Cog thread in a fan-shaped manner into the SMAS layer, from 1 cm anterior to the tragus as the entrance point. The cannula tip should be proceeded 1cm medial to the marionette line.





# Procedures : N-Cog





# Procedures : N-Cog



## LOWER FACE 17

### Submandibular line : N-Cog Multi / N-Cog Spiral

**Entrance point :** ● Lower border of the mastoid process

**End point :** Vertical extension line from marionette line

**Layer :** Subdermal or just above platysma muscle




N-Cog Multi 21GX100 mm

or




N-Cog Spiral 21GX100 mm

To improve submandibular lines, insert N-Cog into the subdermal layer from the mastoid process lower border following the mandibular border. The cannula tip should approach the vertical extension line from marionette line. Insert the 2<sup>nd</sup> and 3<sup>rd</sup> threads below, with a 1 cm distance between each other in a fan-shaped manner.



# Procedures : N-Cog




- Starting point the mastoid process (lower boarder)






# Procedures : N-Cog



full face (double chin)



# Procedures : N-Cog



# Procedures : Y-Ko




- Stimulates patient's intrinsic healing process in order to re-shape the nose.
- Re-Shape your nose in less than 15 minutes and without surgery. Lift the tip of the nose, alignment of the nasal dorsum, and reduction of the nasal wings

## N-Cog Y-KO



Half-Cut cannula (mm)		Thread(mm)	
Gauge	Length	USP	Length
19G	40	0	40
21G	60	2-0	90

30 threads per box. (6 pouches with 5 units in individual sachet)





# Procedures : Y-Ko



## MID FACE 12-1

### Nasal dorsum augmentation, Nasal tip elevation and Nostril reduction : Y-KO

**Entrance point :** ● 0.5 cm below the tip of the nose

**End point :** **Nasal dorsum**\_Glabella, **Columella**\_Anterior nasal spine

**Layer :** **Nasal dorsum**\_Supraperiosteal layer,


**Columella**\_Intercrural space



Y-Ko 21GX60 mm

Y-Ko 19GX40 mm

When performing a nasal dorsum augmentation, insert Y-KO thread through the supraperiosteal layer from 0.5 cm below the nasal tip. Withdraw the cannula once the tip of cannula touches the bone of the glabellar area. For nasal tip elevation and nostril reduction, insert Y-KO thread through the intercrural space of columella from 0.5 cm below the tip of the nose. Withdraw the cannula once the cannula tip touches the anterior nasal spine.




# Procedures : Y-Ko





NEEDLE Gauge: 21G NEEDLE Length: 60mm USP 2-0

Nasal ridge lift



NEEDLE Gauge: 19G NEEDLE Length: 40mm USP 0-0

Nasal tip treatment



1. Nasal ridge insertion : **Insert N-Cog Y-Ko 21G 60mm** needle between 2 Alar cartilage, thought Lateral cartilage to Nasal bone. When you follow the bone, you will meet resistance force before Nasal bone. Upraise the needle for smooth insertion.
2. Nasal tip insertion : **Insert N-Cog Y-Ko 19G 40mm** Use the same insertion point with Nasal ridge insertion, Insert vertically thought the cloume of the nose till you meet the bone.



# Procedures : Y-Ko



# Procedures : Y-Ko



# Procedures : Y-Ko





# Procedures : Derma Spring

+ 더마스프링  
**Derma  
Spring**







# Comparison: type of a thread

	Contour thread	Non-Contour thread
Material	Polypropylene, PDO, PLLA	PDO, PLLA, PCL
Example	Aptos, Silhouette	N-Cog Lift
Action Mechanism	Artificial Lifting Lift up the drooped tissue	Natural Lifting Fix and prevent the drooping of tissue Using thread tissue interaction
Type of thread	Bidirectional cog - cutting, molding, cone	Multiple bidirectional cog - cutting, molding, punching Braided, Weaving Monofilament Spring shape
Complications	Asymmetry Migration and protrusion Early drooping	Edema Purpura




# Comparison: Material of a thread

Material	Polydioxanone / PDO	Polylactide / PLA	Polycaprolactone/PCL	Polypropylene / PP
Molecular structure	 <p><i>p-dioxanone</i> → <i>polydioxanone (PDO)</i></p>			
Absorption	Absorbable	Absorbable	Absorbable	Non-absorbable
Intermediate product	Glyoxylate	Lactic Acid+Hydroxy Acid Acidic substances may change microenvironment of tissue	**6-Hydroxycaproic Acid Acidic substances may change microenvironment of tissue	N/A
End product	CO <sub>2</sub> +H <sub>2</sub> O	CO <sub>2</sub> +H <sub>2</sub> O	CO <sub>2</sub> +H <sub>2</sub> O	Polypropylene
Duration after implantation	6-8 months	1 year	2 year	N/A
Tensile strength	48.3 MPa*	45-70 MPa	23 MPa	40 MPa
Elongation at break	500-600%	85%-105%	4700%	100%
Physical property	Easy to handle	Rigid and brittleness	Flexible	Swollen, after absorbing water
Color	Violet or colorless	Colorless	Colorless	Colorless
Main usage	Suture (Pediatric and Ophthalmic Surgery)	Hard implant (Bone screw)	Drug releasing implant (Contraceptive)	Suture, container (Hernia, Pelvic organ prolapse repair/Syringe)

\* A.K. Sugih et al. / European Polymer Journal 45 (2009) 155–164, Experimental studies on the ring opening polymerization of p-dioxanone using an Al(OiPr)<sub>3</sub>-monosaccharide Initiator system

\*\* Aleksandra et al. / Polish Journal of Chemical Technology, 19(2017), 1, 120—126, Biodegradation of poly(ε-caprolactone) in natural water environments

# Comparison: why PDO?




- Tissue reaction decreased gradually even though implanted material remains\*

	Infiltration		Infiltration content		Angiomatosis		Degree of collagenization	
	AS	EV	AS	EV	AS	EV	AS	EV
14th day	D++	D++	L++ E+ Mp++	L++ E+ Mp++	IMGT ++	C+	CCB++ TCB	CCB+++ TCB
30th day	F++	F+	L+ E+ Mp+	L++ E+ Mp+	C+	C+	CCB++ CoCB	CCB++ TCB
90th day	F+	—	L+	E+	AV+	AV+	CCB++ CoCB	CCB+ TCB

AV, arterioles and venules; C, capillaries; CCB, cellularity of collagen bundles; CoCB, coarse collagen bundles; D, diffuse; E, eosinophils; F, focal; IMGT, initial manifestations of granulation tissue; L, lymphocytes; Mp, macrophages; TCB, thin collagen bundles.

- Vulnerable cogs against force



\* G. Sulamanidze et al. / Journal of Cosmetic Dermatology 17(2018) 731–735, Pathomorphological criteria of use efficiency of resorbable and permanent implant in aesthetic medicine and cosmetic dermatology.

# Comparison: why PDO?

## Physical lifting effects

- Long Duration means good for lifting?
- **PDO** < PLA < PCL < PP
- Long duration does not guarantee the best result.
- After some time passed from thread insertion, physical lifting effect will gone after **6~12 months** for all type of the thread.
- PP is non-absorbable thread, looks have a good physical effect from a start, but it will dropped down in the same time of PDO thread.
- Stay long and make long **inflammantory response** might **interrupt general function** of the subject tissues and might cause **granuloma** with harsh foreign body reactions.

Research

JAMA Facial Plastic Surgery | Original Investigation

### Assessment of Scaffolding Properties for Chondrogenic Differentiation of Adipose-Derived Mesenchymal Stem Cells in Nasal Reconstruction

Serban San-Marina, MD, PhD, Ayushman Sharma, MD, Stephen G. Voss, MS, Jeffrey R. Janus, MD, Grant S. Hamilton III, MD

[Supplemental content](#)

**IMPORTANCE** Nasal reconstruction in patients who are missing a significant amount of structural nasal support remains a difficult challenge. One challenge is the deficiency of cartilage left within the nose as a consequence of rhinectomy or a midline destructive disease. Historically, the standard donor source for large quantities of native cartilage has been costal cartilage.

**OBJECTIVE** To enable the development of protocols for new mesenchymal stem cell technologies as alternative procedures with reduced donor site morbidity, risk of infection and extrusion.

**DESIGN, SETTING, AND MATERIALS** We examined 6 popular scaffold materials in current practice in terms of their biodegradability in tissue culture, effect on adipose-derived mesenchymal stem cell growth, and chondrogenic fate commitment. Various biomaterials of matching size, porosity, and fiber alignment were synthesized by electrospinning and overlaid with rabbit adipose-derived mesenchymal cells in media supplemented or not with chondrogenic factors. Experiments were performed in vitro using as end points biomarkers for cell growth and chondrogenic differentiation. Polydioxanone (PDO), poly-3-hydroxybutyrate-co-3-hydroxyvalerate (PHBV), PHBV-polycaprolactone, poly(L-lactide-co-caprolactone), poly(lactic-co-glycolic acid), and polystyrene scaffolds of 60% to 70% porosity and random fiber alignment were coated with poly(L)-lysine/laminin to promote cell adhesion and incubated for 28 days with  $2.5 \times 10^5$  rabbit adipose mesenchymal cells.

**MAIN OUTCOMES AND MEASURES** Cell growth was measured by fluorometric DNA quantitation and chondrogenic differentiation of stem cells by spectrophotometric sulfated glycosaminoglycan (sGAG) assay. Microscopic visualization of cell growth and matrix deposition on formalin-fixed, paraffin-embedded tissue sections was performed, respectively, with nuclear fast red and Alcian blue.

**RESULTS** Of 6 scaffold materials tested using rabbit adipose mesenchymal cells, uncoated scaffolds promoted limited cell adhesion but coating with poly(L)-lysine/laminin enabled efficient cell saturation of scaffold surfaces, albeit with limited involvement of scaffold interiors. Similar growth rates were observed under these conditions, based on DNA content analysis. However, PDO and PHBV/PCL scaffolds supported chondrogenic fate commitment better than other materials, based on soluble sGAG analysis and microscopic observation of chondrogenic matrix deposition. The mean (SD) sGAG scaffold values expressed as fold increase over control were PDO, 2.26 (0.88), PHBV/PCL, 2.09 (0.83), PLCL, 1.36 (0.39), PLGA, 1.34 (0.77), PHBV, 1.07 (0.31), and PS, 0.38 (0.14).

**CONCLUSIONS AND RELEVANCE** These results establish materials, reagents, and protocols for tissue engineering for nasal reconstruction using single-layer, chondrogenically differentiated, adipose-derived mesenchymal stem cells. Stackable, scaffold-supported, multisheet bioengineered tissue may be generated using these protocols.

**LEVEL OF EVIDENCE** NA.

JAMA Facial Plast Surg. doi:10.1001/jamafacial.2016.1200  
Published online October 13, 2016.

**Author Affiliations:** Department of Otolaryngology, Mayo Clinic, Rochester, Minnesota.


**Corresponding Author:** Serban San-Marina, MD, PhD, Department of Otolaryngology, Mayo Clinic, 200 1st St SW, Rochester, MN 55902 (sanmarina.serban@mayo.edu).

Copyright 2016 American Medical Association. All rights reserved.

Downloaded From: <http://jamanetwork.com/> by See Choi on 11/28/2016

# Comparison: why PDO?


- Which one makes more Collagen?
- PLA < **PDO**
- As we do material experiment for PDO and PLA, PDO makes **7 times more collagen** than PLA.




**FIGURE 3** Changes in expression of TGFβ1, TGFβ2, TGFβ3, collagen type 1, and collagen type 3 by real-time PCR at 1, 2, and 12 wk after injection

KIM ET AL.

JCD  
Journal of  
Cutaneous Medicine and Surgery  
WILEY 13



**FIGURE 1** Histologic evaluation of polydioxanone (PDO). A, PDO is well noted in a 2-wk specimen. B, PDO is still noted with surrounding macrophages and histiocytes in a 12-wk specimen. PDO material is torn down into pieces by fibrous tissues. (A-B; hematoxylin and eosin x100)



**FIGURE 2** Histologic evaluation of poly-L-lactic acid (PLLA). A, PLLA is well noted in a 2-wk specimen. B, PLLA is still noted with a number of macrophages and multinucleated giant cells in a 12-wk specimen. (A-B; hematoxylin and eosin x100)

### 3.3 | Tissue growth factor and collagen changes

Analysis of molecular changes revealed that Col1α1, Col3α1, TGF-β1, TGF-β2, and TGF-β3 showed increases 2 weeks after injection but showed decreases 12 weeks after injection with both PDO and PLLA. Difference value between 1 and 2 week increased more in PDO than PLLA in case of Col1α1, Col3α1, TGF-β2, and TGF-β3. However, TGF-β1 showed more increase in PLLA than PDO (Figure 3, Table 1).

### 4 | DISCUSSION

Currently, fillers for facial volume augmentation are largely divided into two types, hyaluronic acid (HA) and biostimulatory fillers


including CaHA, PCL, and PLLA. Because HA draws water, it has an instant volume-filling effect but a minimal biostimulatory effect; in contrast, PLLA has no instant filling effect but a continuous biostimulatory effect for at least 2 years.<sup>8</sup> Recently, one comparative study suggested that these biostimulatory dermal fillers, such as PCL, are more effective than HA filler for nasolabial fold treatment. Therefore, many high molecular weight polymer-containing fillers that have biostimulatory effects appear continuously today.<sup>7</sup>

Poly-L-lactic acid has biostimulatory effects through macrophages, (myo-)fibroblasts, and collagen that orchestrate capsule formation. The reaction starts with subclinical inflammation that continues to microparticle encapsulation and then to collagen production by the host.<sup>20</sup> PLLA itself degrades into carbon dioxide and water, and the surrounding collagen production continues even after



# Comparison: why PDO?

- **Chemical feature Safty**
- PP < PCL < PLA < **PDO**
- PLA leaves **lactic acid and hydroxy acid** as a intermediate product during biodegradation.
- This might cause hardship to get the correct body test result such as CPR test of subject during thread remains.
- Meanwhile, **PDO will completely absorbed** and excreted **mainly in urine**.



# Conclusion : why PDO?



**PDO** is the most safe and effect material for lifting.



# Remarks

- **High quality PDO**
- **100% biocompatible material**
- **Fast and painless procedures**
- **No General anesthetic needed**
- **Minimal invasive procedure with the most natural results**

# Thank you

**N-FINDERS** | #1211, 55, Digital-ro 33-gil, Guro-gu, Seoul, Korea [08376]  
Tel : +82-70-4801-0558 E-mail : [nfinders@nfinders.com](mailto:nfinders@nfinders.com) [www.nfinders.com](http://www.nfinders.com)