

NEW

Ultra Type

MONALISA

Hyaluronic Acid Body Filler

- Attractive volume to your body
- Proven long-term safety and stability
- Highly pure hyaluronic acid

*Long-term
Stability*




GENOSS
For Patients & Doctors

Product Specification

HA Concentration	24 mg/mL
Lidocaine Content	0.3 %
Volume	10.0 mL
Particle Size	Ave. 900 μ m
Endotoxin	< 0.5 EU/mL


Safety¹

Minimize residual amounts of Endotoxin and BDDE (cross-linker) with optimized purification process.



* 1 GENOSS Data on file

Comparison of Key Properties with Competitors² Higher elasticity and longer duration



* 2 GENOSS Data on file

Volumizing Capacity Marked difference

